

Baxter

SYNCLARA COUGH SYSTEM



THE MOST COMPLETE
COUGH SOLUTION



GOING THE EXTRA MILE FOR PATIENT HEALTH AND COMFORT

For patients with a weak, ineffective cough, cough simulation can be a lifesaving therapy.

At Baxter, we set out to make that therapy as effective and comfortable as possible. The **Synclara** Cough System is the result.

The **Synclara** Cough System uses Mechanical Insufflation-Exsufflation (MIE) technology to clear secretions from the upper airways.


We've added features that enhance individualization of therapy and patient comfort, such as Patient Synchrony and the novel Positive Airway Pressure (PAP) on Pause. These patient-centric features, plus the support service Baxter is known for, make the **Synclara** Cough System the most complete cough solution—for use from hospital to home.

**WE BELIEVE YOUR PATIENTS
DESERVE NOTHING LESS.**

Advanced Patient Synchrony sensing technology automatically adjusts the cough cycle to patients' breathing patterns.

INDIVIDUALIZED CARE


Baxter designed the **Synclara** Cough System to be easily individualized, so you can optimize treatment for a broad range of patients with complex and changing needs, in acute and homecare settings.



DELIVERY OPTIONS


Choose from automatic or manual therapy delivery modes.

NOTE: Manual mode use in home setting determined by Hillrom certified trainer at time of training.




PROGRAMMABLE CYCLES

Up to 10 programmable cycles allow for individualized therapies for your patients.



PAP ON PAUSE

Unique Positive Airway Pressure on Pause feature delivers air between cough cycles, keeping the lungs open and helping patient comfort.



FLUTTER FEATURE

Oscillations delivered to the upper airways during the cough cycle and pause phases help move mucus.

FLEXIBLE DESIGN

Adherence to treatment regimens is critical to improving outcomes, so the **Synclara** Cough System was designed to make therapy as easy and convenient as possible for patients and caregivers.







- 1 PORTABLE**
Small and light, adaptable to virtually any setting
- 2 INTUITIVE TOUCHSCREEN**
Simple to navigate and operate, with preprogrammed protocols ready at the push of a button
NOTE: Manual mode use in home setting determined by Baxter certified trainer at time of patient training
- 3 CIRCUIT WITH MULTIPLE INTERFACES**
Use with mask, mouthpiece or tracheostomy
- 4 ACUTE CARE STAND**
Space-efficient, and easy to move throughout the hospital
- 5 SMART-FILTER WITH RFID**
The **Synclara** Cough System won't run without sensing the **Smart-Filter**
- 6 BATTERY OPERATED (OPTIONAL)**
Take patients' therapy anywhere, with rechargeable battery or direct wall operation

CARE CONNEX PROGRAM FOR SUPPORT



Baxter created the **CARE Connex** program to deliver comprehensive services at every step of the patient journey, offering new efficiencies for your practice and helping improve the therapy experience for your patients. We work with over a thousand leading local, regional and national companies to include commercial health plans and networks, Medicaid, Medicare, TRICARE and the Department of Veterans Affairs.

-  Insurance provider services
-  24/7 clinical care with multilingual customer service and product support
-  600+ in-home licensed clinicians
-  Continuing education courses
-  Clinical in-home assessments, trainings, outcomes monitoring and reporting
-  E-Prescribing available (Parachute, Epic, Cerner, Allscripts and others)

FOR FURTHER INFORMATION ABOUT PRODUCTS AND SERVICES

Please contact your Baxter sales representative or
call Baxter customer service at **1-800-426-4224**.

Rx Only. For safe and proper use of product mentioned herein, please refer to the Instructions for Use or Operator Manual.

Baxter.com

Baxter International Inc.
One Baxter Parkway / Deerfield, Illinois 60015

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Baxter, CARE Connex, Hillrom, Smart-Filter and Synclara are trademarks of Baxter International Inc. or its subsidiaries.

US-FLC172-220026 (v3.0) 02/2024