

TECHNICAL SPECIFICATIONS

Bed	
Low Position - Floor to Center of Deck	13" (33 cm)
Low Position - Floor to Edge of Deck	17" (43.2 cm)
High Position - Floor to Center of Deck	30" (76.2 cm)
High Position - Floor to Edge of Deck	34" (86.4 cm)
Maintains low height when placed in chair	Yes
Headboard/push handles maintain same height regardless of patient surface height	Yes
Sleep Surfaces	
Width	35.5" (90.2 cm)
Length - fully extended	86" (218 cm)
Length - fully retracted	75" (191 cm)
Radiolucent Sleep Deck	No
Surface Weight Limit	500 lbs
Overall Width	
Siderails Stowed	37" (94.0 cm)
Siderails Up	40" (102.0 cm)
Overall Length	
Fully Extended	94.5" (240.0 cm)
Fully Retracted	82.5" (210.0 cm)
Bed Angles	
Head Section	0°to 65°
Knee Section	0° to 16°
Foot Section	0° to -27°
Max Trend/Reverse Trend	15° / 10°
Head and Trend Angle Indicators	Yes
Caster & Braking System	
Caster Diameter	5" (12.7 cm)
Four Wheel Braking	Yes
Brakes Lock Both Rolling and Swiveling	Yes
Scale System	
Accuracy	1% of patient weight
Electrical Characteristics	
Backup Battery Power	Standard
Standards Test to	UL2601 & IEC 60601-1

STANDARD FEATURES

- Scale System
- Side Patient Bed Controls
- SafeView® Alerts*
- Bed Exit Alarm System

* Based on rental availability



For more information, please contact your Hillrom sales representative at 800-638-2546 or order online at rental.hillrom.com.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. 170419 rev 3 28-APR-2021 ENG – US



VersaCare® Bed
Rental offering

SIMPLE, YET SOPHISTICATED

The VersaCare® bed is designed to enable excellent care, especially for mid- to high-acuity patients. We've integrated a host of safety-promoting, leading-edge technologies into one versatile healing system — all controlled through a simple, intuitive interface.



COMPREHENSIVE SAFETY-PROMOTING MEASURES

The VersaCare bed is designed to be able to adjust height, including low positions, to help make it easy for patients to get in and out of bed. Its LowChair® position allows a patient to sit upright while keeping the bed close to the floor.

SafeView® Alerts offer caregivers a visual indication of the safety condition of the bed, including the status of siderails, brakes, HOB and bed height.

If your patient's head-of-bed is lowered below 30 degrees, the Head-of-Bed alarm alerts you of the need to respond.

The 3-Level Bed Exit Alarm with advanced sensitivity detects whether the patient is lying in bed, sitting on the edge of bed, or is out of bed.

- 1 BOOST® ONE-BUTTON FEATURE**
Helps patients move toward the head of the bed
- 2 POINT-OF-CARE® CONTROL PANEL**
Uses easy to-read digital displays and simple icons
- 3 FULLY-FEATURED LOW BED**
Adjusted to the optimal height for patient care with full functionality in the lowest position
- 4 HANDSFREE® FOOT PEDAL SYSTEM**
Allows caregivers to raise and lower the head of the bed and the height of the bed without using their hands
- 5 LOWCHAIR® POSITION**
Allows your patient to be placed in an upright sitting position while keeping the bed low
- 6 FLEXAFOOT™ RETRACTABLE FOOT MECHANISM**
Allows the caregiver to customize the length of the bed up to 11"