



Welch Allyn
Trimline™ Blood Pressure Cuffs



Welch Allyn Trimline™ Blood Pressure Cuffs

Standardize on a proven, dependable cuff technology solution

Why Choose Trimline?

Trimline technology offers reliable cuff connection options that work with virtually any blood pressure device. So it's easy for hospitals to implement a single-patient-use cuff model—and help reduce the risk of cuff-related cross-contamination—while having no impact on everyday clinical workflows.



How can Trimline help you?



Infection Control

Enhance Patient Safety

Blood pressure cuffs are known carriers of pathogens like clostridium difficile.* Hospitals that standardize on Trimline disposable cuffs in a single-patient-use model can reduce the risk of cuff-related cross-contamination, and enhance patient safety in the process.



Clinicians

Build Staff Satisfaction

Using Trimline cuffs, clinicians familiar with traditional tube-and-connector configurations will experience minimal change with how they work today—while your hospital reaps big benefits of enhanced patient safety and improved financial outcomes.



Purchasing

Improve Financial Outcomes

Standardizing on a single blood pressure cuff supplier can help hospitals reduce SKUs, streamline cuff inventory, and save the associated costs.



Biomed Engineers

Enhance Patient Safety

Trimline cuffs use a proven, dependable technology you can count on to work with most blood pressure devices you have in-house.

* Implementation of Disposable Blood Pressure Cuffs as a Novel Approach to Reduce Fomite Transmission of Healthcare-Associated (HCA) Clostridium difficile Infection (CDI) in a Community Hospital or Twice Implemented is Once Credible, American Journal of Infection Control, June 2009.



Welch Allyn Trimline™ Disposable BLOOD PRESSURE CUFFS

- Traditional disposable cuff option that can help reduce the risk of cross-contamination when used in a single-patient-use model
- Withstands extended-use applications
- Available in 'Isolation' color for high-acuity areas
- Exterior sizing guide and balanced bladder design promote accurate sizing—and help to ensure accurate readings
- Full range of sizes available, meets all applicable AAMI standards*



Welch Allyn Trimline™ Reusable BLOOD PRESSURE CUFFS

- Traditional multi-use reusable cuff option that can be cleaned between uses and patients
- Can be used in conjunction with Trimline Disposable Blood Pressure Cuffs to cover just about all hospital settings
- Color coding, exterior sizing guide and balanced bladder design promote accurate sizing—and help to ensure accurate readings
- Two-year warranty



Why standardize with Welch Allyn Trimline™ Cuffs?

- Enables your facility to implement a single-patient-use workflow to help reduce the risk of cuff-related cross-contamination
- Typical hospital savings of 10%-15% on the cost of cuffs
- Helps reduce monitor downtime due to missing specialty cuffs
- Frees up nursing, purchasing and biomedical engineering staff from having to search for cuffs, match up cables and connectors, place emergency orders and manage multiple brands



		Soft Cuffs											
		One-Tube					Two-Tube				Two-Tube with Bulb/Valve		
Size (20 cuffs/case)	Range	Bare Tube	Screw (threaded)	Male Bayonet Rectus	Female Subminiature	Male Subminiature	Bare Tubes	Screw (threaded)	Female Subminiature	Mated Subminiatures	Bare Tube and Clamp	Screw and Clamp	Pocket Aneroid
Infant	9.0-14.8 cm	39239	39090	39045	39026	39127	39008	39010	39163	39030	39171	39119	39241
Small Child	11.8-18.0 cm	39117	39135	39024	39058	39020	39165	39128	39046	39151	39188	39235	39262
Child	13.8-21.5 cm	39083	39091	39100	39007	39044	39154	39011	39249	39031	39115	39120	39186
Small Adult	20.5-28.5 cm	39027	39092	39085	39238	39242	39107	39012	39213	39032	39050	39121	39006
Adult	27.5-36.5 cm	39004	39073	39048	39080	39105	39246	39023	39268	39033	39167	39122	39130
Adult Long	27.5-36.5 cm	39042	39093	39019	39248	39270	39028	39053	39081	39034	39126	39172	39256
Large Adult	35.5-46.0 cm	39124	39074	39230	39039	39261	39103	39014	39149	39035	39254	39021	39185
Large Adult Long	35.5-46.0 cm	39240	39075	39060	39237	39158	39047	39071	39018	39036	39125	39229	39166
Thigh	45.0-56.5 cm	39001	39076	39015	39236	39102	39148	39072	39156	39037	39029	39131	39123
Common Monitor Application		-	Welch Allyn Spot™ Vital Signs, Atlas	Philips, Spacelabs, Datascopie, Siemens, Colin, Prodigy II	Welch Allyn Protocol, Propaq*	Nihon Kohden	-	GE/Critikon Dinamap	GE Marquette	GE Datex-Ohmeda	Manual	Manual	-

		Soft Isolation Cuffs				Vinyl Cuffs						
		Connector type				One-Tube		Two-Tube			Two-Tube with Bulb/Valve	
Size (20 cuffs/case)	Range	Male Bayonet Rectus	Screw (threaded) with Clamp	Bare Tube and Bulb with Clamp	Pocket Aneroid and Bulb	Screw (threaded)	Male Bayonet Rectus	Bare Tubes	Screw (threaded)	Female Subminiatures	Screw and Clamp	Pocket Aneroid
Infant	9.0-14.8 cm	39226	39206	-	-	99851	99866	99842	99850	99629	-	-
Child	13.8-21.5 cm	39225	39205	-	-	99853	99867	99816	99852	99627	-	-
Small Adult	20.5-28.5 cm	39224	39204	-	-	99855	99868	99820	99854	99625	-	-
Adult	27.5-36.5 cm	39223	39203	39209	39210	99859	99869	99825	99856	99624	99882	99888
Adult Long	27.5-36.5 cm	39228	39208	39211	-	99858	99870	99811	99857	99623	99886	99887
Large Adult	35.5-46.0 cm	39222	39202	39212	-	99863	99871	99837	99860	99622	99885	99889
Large Adult Long	35.5-46.0 cm	39227	39207	-	-	99862	99872	99804	99861	99621	-	-
Thigh	45.0-56.5 cm	39221	39201	-	-	99865	99874	99833	99864	99620	-	-

Welch Allyn Isolation Stethoscopes (1/pack) | Yellow: 59803 Red: 59804

		Soft Cuff Connector Type					Vinyl Cuff Connector Type					
Size (40 cuffs/case)	Range	One-Tube Female NeoQuik	One-Tube Male Slip Luer	Two-Tube Female NeoQuik	Two-Tube Male Slip Luers	Two-Tube Male Slip Luer and Clamp	One-Tube Female NeoQuik	One-Tube Male Slip Luer	One-Tube Male Micro Luer	Two-Tube Female NeoQuik	Two-Tube Male Slip Luers	Two-Tube Male Micro Luers
#1 Small	3.4-6.3 cm	NEO-1-ITLS	39190	NEO-1-2TLS	39110	39140	NEO-1-ITLV	99801	99789	NEO-1-2TLV	99760	99832
#2 Medium	4.6-8.5 cm	NEO-2-ITLS	39191	NEO-2-2TLS	39111	39141	NEO-2-ITLV	99890	99844	NEO-2-2TLV	99818	99721
#3 Large	5.8-10.8 cm	NEO-3-ITLS	39192	NEO-3-2TLS	39112	39142	NEO-3-ITLV	99729	99741	NEO-3-2TLV	99873	99806
#4 Extra Large	7.0-13.0 cm	NEO-4-ITLS	39193	NEO-4-2TLS	39113	39143	NEO-4-ITLV	99848	99836	NEO-4-2TLV	99875	99714
#5 Infant	8.2-15.3 cm	NEO-5-ITLS	39194	NEO-5-2TLS	39114	39144	NEO-5-ITLV	99750	99707	NEO-5-2TLV	99722	99824

		One-Tube				Two-Tube			Two-Tube with Inflation System
Size (5 cuffs/bag)	Range	Screw	Bayonet	Female Subminiature	Male Subminiature	Screw	Female Subminiatures	Mated Subminiatures	Bulb/Valve and TP connector
Infant (7)	9.0-14.8 cm	3606SCK	3606SHP	3606SMQ	36706	3606CK	3606MQ	3606DX	3306
Child (9)	13.8-21.5 cm	3605SCK	3605SHP	3605SMQ	36705	3605CK	3605MQ	3605DX	3305
Small Adult (10)	20.5-28.5 cm	3604SCK	3604SHP	3604SMQ	36704	3604CK	3604MQ	3604DX	3304
Adult (11)	27.5-36.5 cm	3603SCK	3603SHP	3603SMQ	36703	3603CK	3603MQ	3603DX	3303
Adult Long (11L)	27.5-36.5 cm	3603XLCK	3603XLSHP	3603XLMSQ	36708	3603XLCK	3603XLMQ	3603LDX	3303XL
Large Adult (12)	35.5-46 cm	3602SCK	3602SHP	3602SMQ	36702	3602CK	3602MQ	3602DX	3302
Large Adult Long (12L)	35.5-46 cm	3602XLCK	3602XLSHP	3602XLMSQ	36709	3602XLCK	3602XLMQ	3602LDX	3302XL
Thigh (13)	45-56.5 cm	3600SCK	3600SHP	3600SMQ	36700	3600CK	3600MQ	3600DX	3300

For more information, contact your local Welch Allyn representative or visit us online at www.welchallyn.com/cuffs.

Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361